

Thank you for downloading [Adaptive Use Musical Instruments](#) (AUMI).

After opening AUMI:

- click on the “click here to set the sound folder” and choose the sound preset folder that is included in your AUMI download.
- Cancel the 3 'save' windows.
- Click the box next to ON&OFF in lower yellow box
- Choose a mode (after choosing a mode the video will display in the black box)
- Change the setting in the blue box on the bottom right to your preference.
- For all modes excluding the Relative Movement Mode click on the part of the body you would like to track. A virtual red marker will appear. (Relative Movement Mode does not need the virtual red marker)
- Now you are ready to make music.

Trouble Shooting:

- If using an external camera be sure the camera is plugged in and recognized by the computer before opening AUMI.
- If there is no sound in percussion modes be sure the 'sound preset' folder is chosen. If already using the patch you can retrieve the folder by clicking in the blue box next to the Mode Name on 'click here to change sound folder'
- If there are further questions please send an email to the AUMI Technical Assistant at jackie.heyen@deeplisting.org. Jackie is also available for skype troubleshooting.

This is beta software and we are constantly working on updates. We will send a message to the email you supplied when there is an update available for download.

Stay connected:

[Donate](#)

[AUMI Website](#)

[Facebook](#)